



PRODUCT INFORMATION

PRODUCT: TRANS-MEDIC

PART NO.: M3612C

FUNCTION AND APPLICATION: Compounded to reduce friction in automatic transmissions and to replenish depleted ATF additives.

OUTSTANDING PROPERTIES: Trans-Medic eliminates transmission slippage, unnecessary downshifting and excessive gear whine. Helps to smooth rough shifting and sluggish action. Works in all automatic transmissions and mixes with all regular fluids. Unaffected by extreme temperatures. Non-corrosive, non-clogging and non-foaming.

USE DIRECTIONS: Check level of transmission fluid according to the instructions from the owner's manual. If the level is correct, drain off 350 ml, of transmission fluid. Add the contents of the container to the transmission. Recheck fluid level. Drive the car normally. The effects of TRANS-MEDIC® will become apparent during 160 km of regular driving. Extremely troublesome cases may require additional treatment with TRANS-MEDIC®. Especially effective when used in conjunction with transmission fluid/filter changes.

PRODUCT LIMITATIONS: A transmission with mechanical defects, seal breakage or excessive sludge accumulation cannot be corrected with Trans-Medic.

PRECAUTIONS: View MSDS at www.gunk.ca.

STORAGE & HANDLING: Store in a cool place below 50°C. Keep away from acids, oxidizing agents and ignition sources. Use with adequate ventilation. Keep away from children and animals.

DISPOSAL: Safeguard our environment. Dispose of container, residue and ATF drained in accordance with all federal, provincial and local laws and regulations.

PACKAGING: 350ml plastic bottle, packed 12 per case.

U.P.C. CODE: 0-60959-02873-9

WEIGHT PER CASE: 11 lbs. 5 kg

CASE SIZE: 10 x 9 x 6 1/2