



**PRODUCT INFORMATION**

**PRODUCT: VALVE-MEDIC**

**PART NO.: M3712C**

**FUNCTION AND APPLICATION:** A scientific blend of solvents compounded to free up noisy valve lifters by dissolving varnish and gum deposits resulting from city traffic and short trips.

**OUTSTANDING PROPERTIES:** One or two containers of VALVE-MEDIC® added to crankcase oil will generally free-up sticking valves and/or valve lifters. Product is quick acting, providing answer to troubled vehicle owners quickly. Easy to use. No tools or mechanical knowledge is required.

**USE DIRECTIONS:** Run engine to its normal operating temperature. Cut off engine. Pour VALVE-MEDIC® into crankcase and restart engine. The VALVE-MEDIC® will usually quiet the engine in 5 to 10 minutes of driving. For stubborn cars, use two containers. If trouble persists, mechanical repair is indicated.

**PRODUCT LIMITATIONS:** VALVE-MEDIC® will not fix worn or damaged parts. If oil is neglected or crankcase anti-pollution system ventilation is defective, noise will probably return.

**PRECAUTIONS:** View MSDS at [www.gunk.ca](http://www.gunk.ca).

**STORAGE & HANDLING:** Store in a cool place, below 50°C away from ignition sources, oxidizing agents and acids. Keep away from children and animals.

**DISPOSAL:** Safeguard our environment. Dispose of container, residue and excess product in accordance with all federal, provincial and local laws and regulations.

**PACKAGING:** 350ml plastic bottle with C/R cap, packed 12 per case.

**U.P.C. CODE:** 0-60959-13371-6

**WEIGHT PER CASE:** 10 lbs. 4.54 kg

**CASE SIZE:** 10 x 7 x 9